

# MB系列大量程天平

- 外观时尚简约, LCD 液晶显示屏
- 可在全量程范围内去皮、清零
- 超载, 欠载显示, 故障报警
- 有计数等多种功能可供选择
- 交直流两用
- 标配RS232 接口
- 可选配天平专用打印机
- 7种单位选择



型号参数	最大称量值 (kg)	可读性 (g)	重复性 (g)	线性 (g)	秤盘尺寸 (mm)
MB10K	10	1	±1	±2	285x235
MB20K	20	1	±1	±2	285x235
MB30K	30	1	±1	±2	285x235
MB10K-1	10	0.1	±0.1	±0.2	285x235
MB20K-1	20	0.1	±0.1	±0.2	285x235
MB25K-1	25	0.1	±0.1	±0.2	285x235

